Absorbable Sutures

Atramat® PGC25 Surgical Suture

General Characteristics

PGC25 Atramat® is a monofilament, synthetic absorbable copolymer of glycolide and caprolactone [poly(glycolide-co-epsilon-caprolactone)].

Characteristics

- High tensile strength.
- PGC25 Atramat® is indicated in surgeries for the closure of surgical layers, which requires a suture material with a useful resistance period of approximately 14-18 days.
- Absorption by hydrolysis.
- Monofilament, uniform diameter.
- Violet or colorless.
- Sterilization by Ethylene Oxide.
- Expiration 3 years.

Absorption

- Loss of tensile strength: 50% at 14 days, 0% at 28 days.
- Absorption is produced by progressive hydrolysis and is complete after 110 days.

Benefits

- Plastic knot coupling providing excellent knot stability.
- Perfect adaptation to the edges of the wound without altering its irrigation due to its good flexibility.
- Elastic and very easy to handle.
- Easy recovery from memory effect
- No capillarity.
- Excellent visibility in the surgical field.
- High resistance to:
 - Traction.
 - Knot breakage.
- Adequate tissue support.
- Predictable absorption rate.
- Smooth passage through tissue.

Applications

It is indicated for:

- The soft approach.
- Suturing limbs, anatomical structures with cavities and hollow liquid-filled organs.
- General and digestive surgery.
- Gynecology and obstetrics.
- Urology.
- Plastic surgery.
- Ligatures.

Knotting Recommendations

It is recommended to use thick gauges (2 - 1 - 0 - 2/0) USP. Make a surgeon's knot, plus a simple knot. This allows the tension to be adjusted. Once adjusted, it is advisable to add a square knot.

In fine gauges (3/0 - 4/0 - 5/0 - 6/0) USP, two initial square knots are recommended with two or three additional fixing knots leaving the ends 4 or 5 mm in length, depending on the case and the preference of the surgeon.

High Quality Needles

Excellent quality 300 grade stainless steel needles with an optimum geometric design and a technologically innovative coating, resulting in an ideal needle.

ATRAGLIDE Needles: A combination of geometry, sharpness and technology with outstanding performance facilitating penetration and smooth passage through tissue.

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