



282466



282467



282468

NSK Contra Angle 4:1

Cat. No

Using a contra angle, compared to a straight nose cone, gives the vet the possibility of turning and locking the angle head to reach difficult areas of the mouth.

- Reduces polishing from 40,000 to 10,000 rpm
- To be used with contra angle head
- The contra angle and heads are fully autoclavable at 135°C
- Can be used on all air-driven micromotors
- 1-year warranty

NSK contra angle 4:1

282466

NSK contra angle head for snap-on cups

282467

NSK contra angle head w/swing latch for mandril cups

282468

See page 35 for Straight Nose Cones.



282370



271127



271121



271080



282487

Polishing

Cat. No

The first and only ergonomically designed prophy angle with a 90° reciprocating “back and forth” rotating motion that provides continuous contact and pressure on the surface of the tooth without catching muzzle hair or whiskers.

This patented oscillating motion:

- Eliminates splatter of prophy paste
- Avoids cross-contamination
- Eliminates frictional heat, making it friendlier on the surface of the tooth while protecting enamel

KRUUSE oscillating prophy angle with cup, 100/pk

282370

KRUUSE prophy angle, metal, for straight handpiece

271127

Polishing cups, snap on, blue, 100/pk

271121

Polishing cups, screw-in, blue, 100/pk

282487

Polishing cups latch style, 12/pk

271080



282372



282377



282373



282378

Prophy Paste

Cat. No

Prophy paste free of Xylitol and fluoride. Fast stain removal and polishing. Easy to apply, low splatter, rinses easily. Use coarse grit when heavy stains need to be removed.

Prophy paste, coarse grit, red, mint , 2 g, 100/pk

282372

Prophy paste, medium grit, green, mint , 2 g, 100/pk

282373

Prophy paste, coarse grit, red, mint , 255 g

282377

Prophy paste, medium grit, green, mint , 255 g

282378